

One Step Methadone Urine Test

Catalog No. See Pouch Label

One Step Methadone Urine Test is a rapid one step test for the qualitative detection of Methadone and its principal metabolites in human urine at specified cut-off level. For in vitro diagnostic use only. For professional use only.

INTENDED USE

One Step Methadone Urine Test is a lateral flow chromatographic immunoassay for the detection of Methadone in human urine at the cut-off concentration of 300 ng/ml. This assay provides only a qualitative, preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are used.

SUMMARY

Methadone is a synthetic analgesic drug that is originally used in the treatment of narcotic addicts. Among the psychological effects induced by using methadone are analgesia, sedation and respiratory depression. Overdose of methadone may cause coma or even death. It is administered orally or intravenously and is metabolized in the liver and excreted in urine as methadone, EDDP, EMDA and methadol. The kidneys are a major route of methadone excretion. Methadone has a biological half-life of 15 to 60 hours.

PRINCIPLE

One Step Methadone Urine Test is a competitive immunoassay that is used to screen for the presence of Methadone in urine. It is chromatographic absorbent device in which Methadone and its metabolites in a sample competitively combined to a limited number of anti-methadone monoclonal antibody (mouse) conjugate binding sites.

When add four drops of the urine specimen to the sample well, the urine is absorbed into the device by capillary action, mixes with the methadone monoclonal antibody conjugate, and flows across the pre-coated membrane. When sample drug levels are zero or below the target cut off (the detection sensitivity of the test), anti-methadone monoclonal antibody (mouse) conjugate binds to the methadone-protein (duck egg) conjugate immobilized in the Test Region (T) of the device. This produces a colored Test line that, regardless of its intensity, indicates a negative result.

When sample drug levels are at or above the target cutoff, the free drug in the sample binds to the methadone monoclonal antibody conjugate preventing the methadone monoclonal antibody conjugate from binding to the methadone-protein conjugate immobilized in the Test Region (T) of the device. This prevents the development of a distinct colored band in the test region, indicating a potentially positive result.

To serve as a procedure control, a colored line will appear at the Control Region (C), where the Goat anti mouse IgG polyclonal antibody immobilized in, if the test has been performed properly.

WARNINGS AND PRECAUTIONS

1. This kit is for external use only. Do not swallow.
2. Discard after first use. The test cannot be used more than once.
3. Do not use test kit beyond expiry date.
4. Do not use the kit if the pouch is punctured or not well sealed.
5. Keep out of the reach of children.
6. Do not read after 5 minutes
7. This kit is for in vitro diagnostic use.

CONTENT OF THE KIT

1. 25 tests per kit., one test in one pouch.
2. One pouch containing a test and a desiccant. The only desiccant is for storage purposes only, and is not used in the test procedures.
3. Leaflet with instructions for use.
4. Dropper.

STORAGE AND STABILITY

Store at 4 ~ 30 °C in the sealed pouch up to the expiration date. Keep away from direct sunlight, moisture and heat.

DO NOT FREEZE.


SPECIMEN COLLECTION AND PREPARATION

Collect a urine sample in the supplied urine cup. Urine specimens may be refrigerated (2-8°C) and stored up to forty-eight hours. For longer storage, freeze the samples (-20°C or below). Bring frozen or refrigerated samples to room temperature before testing. Previously frozen or refrigerated samples should be well mixed before analysis. Cloudy specimens should be centrifuged before analysis. Use only clear aliquots for testing.

TEST PROCEDURE

Test must be in room temperature (18°C to30°C)

1. Open the sealed pouch by tearing along the notch. Remove the test device from the pouch.
2. Hold the one side of the device with one hand. Use the other hand to pull out the cap and expose the absorbent end.
3. Immerse the absorbent end into the urine sample about 10 seconds. Make sure that the urine level is not above the "MAX" line printed on the front of the device.
4. Lay the device flat on a clean, dry, non-absorbent surface.
5. Read the result at 5 minutes. **Do not read after 5 minutes.**



INTERPRETATION OF RESULTS

Positive (+)

A rose-pink band is visible in the control region. No color band appears in the appropriate test region. It indicates a positive result for the corresponding drug of that specific test zone.


Negative (-)

A rose-pink band is visible in the control region and the appropriate test region. It indicates that the concentration of the corresponding drug of that specific test zone is below zero or the detection limit of the test.

Invalid

If a color band is not visible in the control region or a color band is only visible in the test region, the test is invalid. Another test should be run to re-evaluate the specimen. Please contact the distributor or the store, where you bought the product, with the lot number.

Note: There is no meaning attributed to line color intensity or width.



QUALITY CONTROL

Users should follow the appropriate federal state, and local guidelines concerning the frequency of assaying external quality control materials.

Though there is an internal procedural control line in the test device of Control region, the use of external controls is strongly recommended as good laboratory testing practice to confirm the test procedure and to verify proper test performance. Positive and negative control should give the expected results. When testing the positive and negative control, the same assay procedure should be adopted.

LIMITATIONS OF PROCEDURE

1. This test has been developed for testing urine samples only. The performance of this test using other specimens has not been substantiated.
2. Adulterated urine samples may produce erroneous results. Strong oxidizing agents such as bleach (hypochlorite) can oxidize drug analytes. If a sample is suspected of being adulterated, obtain a new sample.
3. This test is a qualitative screening assay. It is not designed to determine the quantitative concentration of drugs or the level of intoxication

PERFORMANCE CHARACTERISTICS

Accuracy

Eighty clinical urine specimens were analyzed by GC-MS and by the One Step Methadone Test. Each test was read by three viewers. Samples were divided by concentration into four categories: less than half the cutoff, near cutoff negative, near cutoff positive, and high positive. Results were as follows:

Viewer A:

Result	Less than half the cutoff concentration by GC/MS analysis	Near Cutoff Negative (Between 50% below the cutoff and the cutoff concentration)	Near Cutoff Positive (Between the cutoff and 50% above the cutoff concentration)	High Positive (greater than 50% above the cutoff concentration)
Positive	0	0	17	21
Negative	22	18	2	0

% agreement among positives is 95% (95% Confidence Interval 79.5% - 100%)

% agreement among negatives is 100% (95% Confidence Interval 84.5% - 100%)

Viewer B:

Result	Less than half the cutoff concentration by GC/MS analysis	Near Cutoff Negative (Between 50% below the cutoff and the cutoff concentration)	Near Cutoff Positive (Between the cutoff and 50% above the cutoff concentration)	High Positive (greater than 50% above the cutoff concentration)
Positive	0	0	18	21
Negative	22	18	1	0

% agreement among positives is 97.5% (95% Confidence Interval 82% - 100%)

% agreement among negatives is 100% (95% Confidence Interval 84.5% - 100%)

Viewer C:

Result	Less than half the cutoff concentration by GC/MS analysis	Near Cutoff Negative (Between 50% below the cutoff and the cutoff concentration)	Near Cutoff Positive (Between the cutoff and 50% above the cutoff concentration)	High Positive (greater than 50% above the cutoff concentration)
Positive	0	0	18	21
Negative	22	18	1	0

% agreement among positives is 97.5% (95% Confidence Interval 82% - 100%)

% agreement among negatives is 100% (95% Confidence Interval 84.5% - 100%)

From the results of the above tables, the total results are showed as below:

The average positive agreement is 96.7%

The average negative agreement is 100%

Precision and Sensitivity

To investigate the precision and sensitivity, samples were analyzed at the following concentrations: cutoff - 50%, cutoff - 25%, cutoff, cutoff +25%, and the cutoff + 50%. All concentrations were confirmed with GC-MS. Each concentration was tested using three different lots. Thirty samples were analyzed at each concentration, and each result was read by three viewers, for a total of 90 results per concentration per lot.

Lot 1

Approximate concentration of sample (ng/mL)	Number of determinations	Results Negative/ Positive
150	90	90/0
225	90	75/15
300	90	41/49
375	90	7/83
450	90	0/90

Lot 2

Approximate Concentration of sample (ng/mL)	Number of determinations	Results Negative/ Positive
150	90	90/0
225	90	75/15
300	90	41/49
375	90	7/83
450	90	0/90

Lot 3

Approximate Concentration of sample (ng/mL)	Number of determinations	Results Negative/ Positive
150	90	90/0
225	90	75/15
300	90	41/49
375	90	7/83
450	90	0/90

Specificity and cross reactivity

To test the specificity of the test, the test device was used to test Methadone, drug metabolites and other components that are likely to be present in urine. All the components were added to drug-free normal human urine. These concentrations below also represent the limits of detection for the specified drugs or metabolites.

Component	Concentration (ng/ml)
-----------	-----------------------

Methadone	300 (100%)
Doxylamine	50,000 (0.6%)

Effect of Urinary Specific Gravity

5 urine samples with density ranges (1.000-1.035) are collected and spiked with methadone at 50% below and 50% above cutoff level. One step Methadone urine test was tested in duplicate. The results demonstrate that varying ranges of urinary specific gravity do not affect the test result.

Effect of Urinary PH

The pH of an aliquot negative urine pool is adjusted to a pH range of 4 to 9 in 1 pH unit increments and spiked with methadone at 50% below and 50% above cutoff levels. One step Methadone urine test was tested in duplicate. The result demonstrate that varying ranged of PH do not interfere with the performance of the test.

Interfering substances

Clinical urine samples may contain substances that could potentially interfere with the test. The following compounds were added to drug-free urine, urine with a Methadone concentration 50% below the cutoff, and urine with a Methadone concentration 50% above the cutoff. All potential interferents were added at a concentration of 100 µg/mL. None of the urine samples showed any deviation from the expected results.

Acetaminophen	Loperamide
Acetophenetidin	Mephentermine
N-Acetylprocainamide	Maprotiline
Acetylsalicylic acid	Meperidine
Aminopyrine	Meprobamate
Amityptiline	Methamphetamine
Amobarbital	Methoxyphenamine
Amoxicillin	(±) - 3,4-Methylenedioxyamphetamine
Ampicillin	hydrochloride
L-Ascorbic acid	(±)-3,4-Methylenedioxyamphetamine
DL-Amphetamine sulfate	hydrochloride
Apomorphine	Morphine-3-β-D glucuronide
Aspartame	Morphine Sulfate
Atropine	Nalidixic acid
Benzilic acid	Naloxone
Benzoic acid	Naltrexone
Benzoylcegonine	Naproxen
Benzphetamine	Niacinamide
Bilirubin	Nifedipine
(±) - Brompheniramine	Norcodein
Caffeine	Norethindrone
Cannabidiol	D-Norpropoxyphene
Cannabinol	Noscapine
Chloralhydrate	DL-Octopamine
Chloramphenicol	Oxalic acid
Chlorothiazide	Oxazepam
(±) - Chlorpheniramine	Oxolinic acid
Chlorpromazine	Oxycodone
Chlorquine	Oxymetazoline
Cholesterol	Papaverine
Clomipramine	Penicillin-G
Clonidine	Pentazocine hydrochloride
Cocaethylene	Pentobarbital
Temazepam	Perphenazine
Cocaine hydrochloride	Phencyclidine
Codeine	Phenelzine
Cortisone	Phenobarbital
(-) Cotinine	Phentermine
Creatinine	Trans-2-phenylcyclopropylamine
Deoxycorticosterone	hydrochloride
Dextromethorphan	L-Phenylephrine
Diazepam	β-Phenylethylamine
Diclofenac	Phenylpropanolamine
Diffunisal	Prednisolone
Digoxin	Prednisone
Diphenhydramine	Procaine

Ecgonine hydrochloride	Promazine
Ecgonine methylester	Promethazine
(-) -Ψ -Ephedrine	DL-Propranolol
[1R,2S] (-) Ephedrine	D-Propoxyphene
(L) - Epinephrine	D-Pseudoephedrine
Erythromycin	Quinacrine
β-Estradiol	Quinidine
Estrone-3-sulfate	Quinine
Ethyl-p-aminobenzoate	Ranitidine
Fenoprofen	Salicylic acid
Furosemide	Secobarbital
Gentisic acid	Serotonin
Hemoglobin	Sulfamethazine
Hydralazine	Sulindac
Hydrochlorothiazide	Tetracycline
Hydrocodone	Tetrahydrocortisone, 3-acetate
Hydrocortisone	Tetrahydrocortisone 3- (β-D-glucuronide)
O-Hydroxyhippuric acid	Tetrahydrozoline
p-Hydroxyamphetamine	Thebaine
p-Hydroxymethamphetamine	Thiamine
3-Hydroxytyramine	Thioridazine
Ibuprofen	DL-Tyrosine
Imipramine	Tolbutamide
Iproniazid	Triamterene
(±) - Isoproterenol	Trifluoperazine
Isoxsuprine	Trimethoprim
Ketamine	Trimipramine
Ketoprofen	Tryptamine
Labelalol	DL-Tryptophan
Levorphanol	Tyramine
	Uric acid
	Verapamil
	Zomepirac

BIBLIOGRAPHY OF SUGGESTED READING

Baselt, R.C. Disposition of Toxic Drugs and Chemicals in Man. Biomedical Publications, Davis, CA, 1982.

Ellenhorn, M.J. and Barceloux, D. G Medical Toxicology. Elsevier Science Publishing Company, Inc., New York, 1988

Gilman, A. G., and Goodman, L. S. The Pharmacological Fluids, in Martin WR(ed): Drug Addiction I, New York, Spring – Verlag, 1977.

Harvey, R.A., Champe, P.C. Lippincotts Illustrated Reviews. Pharmacology. 91-95, 1992.

Hawwks RL, CN Chiang. Urine Testing for drugs of Abuse. National Institute for Drug Abuse (NIDA), Research Monography 73, 1986

Hofmann F.E., A Handbook on Drug and Alcohol Abuse: The Biomedical Aspects, New York, Oxford University Press, 1983.

McBay, A. J. Clin. Chem. 33,33B–40B, 1987.

Index of Symbols

	For <i>in vitro</i> diagnostic use only		Store between 4 ~ 30 °C
	Do not reuse		Manufacturer
	Keep away from sunlight		Keep dry