

One Step Oxycodone Urine Test

Catalog No. See Pouch Label

One Step Oxycodone Urine Test is a rapid one step test for the qualitative detection of Oxycodone and its principal metabolites in human urine at specified cut-off level. For in vitro diagnostic use only. For healthcare professional use only.

INTENDED USE

One Step Oxycodone Urine Test is a lateral flow chromatographic immunoassay for the detection of Oxycodone in human urine at the cut-off concentration of 100 ng/ml. This assay provides only a qualitative, preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are used.

SUMMARY

Oxycodone is known as Oxycontin, Roxicodone and is an ingredient of Percodan, Percocet, Roxicet and Tylox. Oxycodone is a semi-synthetic opiate derived from opium. Like other opiates, oxycodone is characterized by its analgesic properties, and the tendency for users to form a physical dependency and develop tolerance with extended use. Oxycodone is usually administered in combination with non-opiate analgesics such as acetaminophen and salicylates for the relief of moderate to severe pain. Oxycodone is a central nervous system depressant that may cause drowsiness, dizziness, lethargy, weakness and confusion. Toxicity in an overdose of oxycodone can lead to stupor, coma, muscle flaccidity, severe respiratory depression, hypotension, and cardiac arrest.

Oxycodone is metabolized by N- and O-demethylation. One of the metabolites, oxymorphone, is a potent narcotic analgesic, while the other, noroxycodone, is relatively inactive. Between 33 to 61% of a single dose of oxycodone is excreted in a 24 hour urine collection and consists of 13-19% free oxycodone, 7-29% glucuronide conjugated oxycodone, 13-14% glucuronide conjugated oxymorphone and an unknown amount of noroxycodone. The detection time window of oxycodone is 1-3 days following use.

PRINCIPLE

One Step Oxycodone Urine Test is a competitive immunoassay that is used to screen for the presence of Oxycodone in urine. It is a chromatographic absorbent device in which Oxycodone and its metabolites in a sample competitively combine to a limited number of antibody-dye conjugate binding sites.

When four drops of the urine specimen to the sample well, the urine is absorbed into the device by capillary action, mixes with the antibody-dye conjugate, and flows across the pre-coated membrane. When sample drug levels are zero or below the target cut-off (the detection sensitivity of the test), antibody-dye conjugate binds to the drug-protein conjugate immobilized in the Test Region (T) of the device. This produces a colored Test line that, regardless of its intensity, indicates a negative result.

When sample drug levels are at or above the target cutoff, the free drug in the sample binds to the antibody-dye conjugate preventing the antibody-dye conjugate from binding to the drug-protein conjugate immobilized in the Test Region (T) of the device. This prevents the development of a distinct colored band in the test region, indicating a potentially positive result.

To serve as a procedure control, a colored line will appear at the Control Region (C), if the test has been performed properly.

PRECAUTIONS

1. This kit is for external use only. Do not swallow.
2. Discard after first use. The test cannot be used more than once.
3. Do not use test kit beyond expiration date.
4. Do not use the kit if the pouch is punctured or not well sealed.
5. Keep out of the reach of children.
6. Do not read after 5 minutes

CONTENT OF THE KIT

1. 25 tests per kit, one test in one pouch.
2. One pouch containing a test, a dropper and a desiccant. The desiccant is for storage purposes only, and is not used in the test procedures.
3. Leaflet with instructions for use.
4. Dropper.

STORAGE AND STABILITY

Store at 4 °C ~30 °C in the sealed pouch up to the expiration date. Keep away from direct sunlight, moisture and heat. DO NOT FREEZE.

SPECIMEN COLLECTION AND PREPARATION


Collect a urine sample in the supplied urine cup. Urine specimens may be refrigerated (2-8°C) and stored up to forty-eight hours. For longer storage, freeze the samples (-20°C or below).

Bring frozen or refrigerated samples to room temperature before testing. Previously frozen or refrigerated samples should be well mixed before analysis. Cloudy specimens should be centrifuged before analysis. Use only clear aliquots for testing.

TEST PROCEDURE

Test must be in room temperature (18°C to 30°C)

1. Open the sealed pouch by tearing along the notch. Remove the test device from the pouch.
2. Hold the one side of the device with one hand. Use the other hand to pull out the cap and expose the absorbent end.
3. Immerse the absorbent end into the urine sample about 10 seconds. Make sure that the urine level is not above the "MAX" line printed on the front of the device.
4. Lay the device flat on a clean, dry, non-absorbent surface.
5. Read the result at 5 minutes. **Do not read after 5 minutes.**



INTERPRETATION OF RESULTS

Positive (+)

A rose-pink band is visible in the control region. No color band appears in the test region. This positive result indicates that the Oxycodone concentration is equal to or higher than the detection limit (100 ng/ml).


Negative (-)

A rose-pink band is visible in the control region and the test region. This negative result indicates that the Oxycodone concentration is zero or below the detection limit (100 ng/ml).

Invalid

If a color band is not visible in the control region or a color band is only visible in the test region, the test is invalid. Another test should be run to re-evaluate the specimen. If test still fails, please contact the distributor or the store, where you bought the product, with the lot number.

Note: There is no meaning attributed to line color intensity or width.



QUALITY CONTROL

Though there is an internal procedural control line in the test device of Control region, the use of external controls is strongly recommended as good laboratory testing practice to confirm the test procedure and to verify proper test performance. Positive and negative control should give the expected results. When testing the positive and negative control, the same assay procedure should be adopted.

LIMITATIONS OF PROCEDURE

1. This test has been developed for testing urine samples only. The performance of this test using other specimens has not been substantiated.
2. Adulterated urine samples may produce erroneous results. Strong oxidizing agents such as bleach (hypochlorite) can oxidize drug analytes. If a sample is suspected of being adulterated, obtain a new sample.
3. This test is a qualitative screening assay. It is not designed to determine the quantitative concentration of drugs or the level of intoxication.

PERFORMANCE CHARACTERISTICS

A. Sensitivity

One Step Oxycodone Urine Test has set the screen cut-off for positive specimens at 100 ng/mL for Oxycodone as a calibrator. The test device has been proved to detect above 100 ng/mL of Oxycodone in urine at 5 minutes.

B. Specificity and cross reactivity

To test the specificity of the test, the test device was used to test Oxycodone, its metabolites and other components of the same class that are likely to be present in urine. All the components were added to drug-free normal human urine. These concentrations below also represent the limits of detection for the specified drugs or metabolites.

Component	Concentration (ng/ml)
Oxycodone	100
Dihydrocodeine	20,000
Codeine	100,000
Hydromorphone	100,000
Morphine	> 100,000
Acetylmorphone	> 100,000
Buprenorphine	> 100,000

Ethylmorphine > 100,000

C. Interfering substances

Considering the complexity of clinical urine specimens and the possibility that various urine specimens contain potentially interfering substances, for example Acetoacetic Aci, Acetone, Albumin etc., we simulated above situations by adding the potentially interfering substances to a certain concentration as specimen. The following components show no cross-reactivity when tested with One Step Oxycodone Urine Test at a concentration of 100 µg/ml.

4-Acetamidophenol	Maprotiline
Acetophenetidin	Meperidine
N-Acetylprocainamide	Meprobamate
Acetylsalicylic acid	Methadone
Aminopyrine	Morphine-3-β-Dglucuronide
Amityryptiline	Nalidixic acid
Amobarbital	Naloxone
Amoxicillin	Naltrexone
BUPicillin	Naproxen
Ascorbic acid	Niacinamide
Apomorphine	Nifedipine
Aspartame	Norcodein
Atropine	Norethindrone
Benzilic acid	D-Norpropoxyphene
Benzoic acid	Noscapine
Benzoylcegonine	D,L-Octopamine
Benzphetamine	Oxalic acid
Bilirubin	Oxazepam
Brompheniramine	Oxolinic acid
Caffeine	Oxymetazoline
Cannabidiol	Papaverine
Cannabinol	Penicillin-G
Chloralhydrate	Pentazocaine
ChlorBUPhenicol	Pentobarbital
Chlordiazepoxide	Perphenazine
Chlorothiazide	Phencyclidine
(±) Chlorpheniramine	Phenelzine
Chlorpromazine	Phendimetrazine
Chlorquine	Phenobarbital
Cholesterol	Phetoin
Clomipramine	L-Phenylephrine
Clonidine	β-Phenylethamine
Cocaine hydrochloride	Phenylpropanolamine
Cortisone	Prednisolone
(-) Cotinine	Prednisone
Creatinine	Procaine
Deoxycorticosterone	Promazine
Dextromethorphan	Promethazine
Diazepam	D,L-Propanolol
Diclofenac	Propiomazine
Diffunisal	D-Propoxyphene
Digoxin	Quinidine
Diphenhydramine	Quinine
Doxylamine	Ranitidine
Ecgonine hydrochloride	Salicylic acid
Ecgonine methylester	Secobarbital
(1R,2S)-(-)-Ephedrine	Serotonin
L-Ephedrine	Sulfamethazine
(-) Y Ephedrine	Sulindac
Erythromycin	Temazepam
β-Estradiol	Tetracycline
Estrone-3-sulfate	Tetrahydrocortisone

Ethyl-p-aminobenzoate	Tetrahydrozoline
Fenfluramine	Δ 9-THC-COOH
Fenoprofen	Thebaine
Furosemide	Thiamine
Gentisic acid	Thioridazine
Hemoglobin	D,L-Thyroxine
Hydralazine	Tolbutamine
Hydrochlorothiazide	Triamterene
Hydrocodone	Trifluoperazine
Hydrocortisone	Trimethoprim
O-Hydroxyhippuric acid	Trimipramine
3-Hydroxytyramine	Tryptamine
Ibuprofen	D, L-Tyrosine
Imipramine	Uric acid
(-) Isoproterenol	Verapamil
Isoxsuprine	Zomepirac
Ketamine	
Ketoprofen	
Labetalol	
Levorphanol	
Loperamide	

BIBLIOGRAPHY OF SUGGESTED READING

1. Baselt, R.C. Disposition of Toxic Drugs and Chemicals in Man. Biomedical Publications, Davis, CA, 1982.
2. Ellenhorn, M.J. and Barceloux, D. G Medical Toxicology. Elsevier Science Publishing Company, Inc., New York, 1988
3. Gilman, A. G., and Goodman, L. S. The Pharmacological Fluids, in Martin WR(ed): Drug Addiction I, New York, Spring – Verlag, 1977.
4. Harvey, R.A., Champe, P.C. Lippincotts Illustrated Reviews. Pharmacology. 91-95, 1992.
5. Hawwks RL, CN Chiang. Urine Testing for drugs of Abuse. National Institute for Drug Abuse (NIDA), Research Monography 73, 1986
6. Hofmann F.E., A Handbook on Drug and Alcohol Abuse: The Biomedical Aspects, New York, Oxford University Press, 1983.
- McBay, A. J. Clin. Chem. 33,33B-40B, 1987.

MEANING OF SYMBOLS ON PACKAGE



Keep away from sunlight



Store between 4°C and 30°C



Keep dry



Do not re-use